尺寸：120mm*105mm

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support.it@taotronics.com(IT)

EC REP

Name: Share Info Consultant Service LLC.

Addr.: Heerdter Lohweg 83, 40549 Düsseldorf.

E-mail: eu-rep@share-info.cn

Tel / Fax: 0049 1767 0057022

NORTH AMERICA

E-mail : support@taotronics.com(US)

support.ca@taotronics.com(CA)

Tel : 1-888-456-8468 (Monday-Friday: 9:00 – 17:00 PST)

US Distributor: SUNVALLEYTEK INTERNATIONAL INC.

Address: 46724 Lakeview Blvd, Fremont, CA 94538

MANUFACTURER

Name: Shenzhen Whakin Innovation Technology Co., Ltd.

Addr.: South of the fourth floor, L bldg., jingtie technology industrial park, no.49 changjiang pu road, heao community, yuanshan street, longgang district, shenzhen city, China.

E-mail: qm@whakin.com

Tel / Fax: 00856 755 2828 3837

Learn more about the EU Declaration of Conformity:

<https://www.taotronics.com/downloads-TT-CL016-CE-Cert.html>


1282

MADE IN CHINA

201909

TT-CL016_V1.4 EU

MODEL: TT-CL016
SAD LIGHT
User Guide

CONTENTS

| | |
|----------|-------|
| English | 01-19 |
| Deutsch | 20-39 |
| Français | 40-59 |
| Español | 60-79 |
| Italiano | 80-99 |

EN

PACKAGE CONTENTS

- ① 1 x TaoTronics SAD Light (Model: TT-CL016)
- ② 1 x Power Adapter
- ③ 1 x User Guide
- ④ 1 x Thank You Card
- ⑤ 1 x Travel Pouch

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

1**EN****PRODUCT DIAGRAM (FIG.1)**

- ① Lamp
- ② Power Button/ Timer Setting
- ③ Light Brightness +
- ④ Light Brightness -
- ⑤ Power Port
- ⑥ Foldable Bracket

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IMPORTANT**Warning**

1. Please keep the lamp away from water or any other liquid to avoid electrical shock.
2. For indoor use only.
3. This appliance is not intended for use by persons (including children) with skin damage, reduced sensory or mental capabilities, or lack of experience and knowledge, unless authorized by a professional physician.
4. Keep away from heat and avoid direct sunlight.
5. Do not perform any repair or maintain while using this product for treatment.
6. Do not try to repair or dismantle the lamp by yourself.
7. Do not drop or knock the lamp on hard surfaces.
8. Portable and mobile RF communication equipment may interfere with medical electrical equipment.
9. Do not place the device in a location where it is difficult to operate the disconnecting device (disconnecting device: Adapter).
10. Using an adapter other than the one provided by the manufacturer may result in failure to work, or even a safety incident.
11. Keep the equipment out of reach of children or infants, avoid strangulation caused by the adapter cable.

Cautions

1. Some people may suffer from light sensitivity, discontinue use immediately if any discomfort occurs.

- Never cover the lamp or place anything on top of it during operation to avoid overheating or fire hazards.
- Consult your health care provider before use of this product if you are under medical care or have suffered from an eye condition before.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components.
- Do not look directly look at the light therapy lamp. It only serves as background light.
- If you have any product problems or report an unexpected situation during use, please unplug the power immediately and call the after-sales service hotline.

⚠ This product contains magnets. Keep a proper distance if you are wearing a pacemaker or any other implants.

INTENDED USE

Your new TaoTronics SAD LIGHT is a remarkable product with many features and benefits:

- Delivers LED light to your home, office, dorm room or apartment.
- Designed with a stable stand for easy placement on a convenient spot.
- Cues your body's own energy enhancers to help you regain focus and concentration.
- UV free LED- the built-in LED light panel delivers the full brightness with no pollution, no radiation and no stroboscopic effects.
- 10,000 lux at 10cm/4" from lamp to lift up your spirit.


SYMBOLS AND MEANINGS

| | | | |
|---|---|---|--|
|  | Serial number of product. |  | The date of manufacture. |
|  | Waste disposal in accordance with WEEE directive requirements. |  | The name and address of the manufacturer. |
|  | Refer to instruction manual / booklet. |  | Start/Off, stand-by button. |
| IP22 | Degree of protection against ingress of water and particulate matter. |  1282 | The device conforms to European Union health, safety, and environmental standards. |
|  | Authorized representative in the European Community. |  | Caution |
|  | Temperature limit |  | Humidity limitation |
|  | Atmospheric pressure limitation |  | Keep away from sunlight |
|  | Keep dry |  | Keep upright during transport |


SPECIFICATIONS

| | |
|---------------------|--|
| Model | TT-CL016 |
| Voltage | AC 100-240V, 50/60Hz DC 12V \equiv 1.0A |
| Power | 12W Max. |
| Storage Temperature | -20°C -50°C / -4°F -122°F |

PREPARING FOR USE

1. Put the lamp on a flat and stable surface using the support foot (Fig.2).

2. Connect to a power source with the included power adapter. The unit will enter standby mode, and the Power Button indicator will light up (Fig.3).

USING YOUR SAD LIGHT

Power On / Off

1. Tap the Power Button once to turn on / off the lamp.
2. Press the + / - Button on the panel to choose the most comfortable brightness for you.
3. If power is not disconnected, the lamp will remember the light brightness of last use when turned on again.

Note: Unplugging the lamp will reset the setting.

Set the Timer

Long press the Power button to set the timer for 30 minutes. The LED Indicator will blink white when the timer is duly enabled. Long press the Power button to disable the timer at any time.

Adjust the Light Brightness

Tap and hold the Light Brightness + / - Button to adjust the light brightness.

TROUBLESHOOTING

| Problem | Possible Cause | Solution |
|-----------------------------------|---|--|
| The lamp does not light up | No power supply | Correctly plug in the power adapter |
| | The lamp is not turned on | Refer to the How to Use section |
| The touch panel is not responsive | Hands are wet | Touch the panel with dry hands |
| | Voltage fluctuation causing panel failure | Disconnect power to restart the lamp |
| LED failure | LED component failed | Please contact the TaoTronics Support Team |

| | | |
|-----------|----------------------|--|
| Bad smell | Power adapter failed | Please contact the TaoTronics Support Team |
| | Short circuit | Please contact the TaoTronics Support Team |

WARRANTY AND SERVICE

In case of any problem, please contact the TaoTronics Customer Care Centre in your country via the contact information listed in this User Guide.

Guidance and manufacturer's declaration - electromagnetic emission - for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission

The model TT-CL016 is intended for use in the electromagnetic environment specified below. The customer or the user of the model TT-CL016 should assure that it is used in such an environment.

| Emissions test | Compliance | Electromagnetic environment - guidance |
|--|------------|---|
| RF emissions CISPR 11 | Group 1 | The Model TT-CL016 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The TT-CL016 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations /flicker emissions IEC 61000-3-3 | Applicable | |

Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity

The Model TT-CL016 are intended for use in the electromagnetic environment specified below. The customer or the user of the Model TT-CL016 should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|--|--|--|---|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrostatic transient / burst IEC 61000-4-4 | ± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines | ± 2 kV for power supply lines 100 kHz repetition frequency | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |

| | | | |
|--|--|--|---|
| Surge IEC 61000-4-5 | ± 0.5 kV, ± 1 kV differential mode line-line | ± 0.5 kV, ± 1 kV differential mode line-line | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <p>0 % UT (100 % dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0 % UT (100 % dip in UT) for 1 cycle at 0°</p> <p>70 % UT (30 % dip in UT) for 25/30 cycles at 0°</p> <p>0 % UT (100 % dip in UT) for 250/300 cycle at 0°</p> | <p>0 % UT (100 % dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0 % UT (100 % dip in UT) for 1 cycle at 0°</p> <p>70 % UT (30 % dip in UT) for 25/30 cycles at 0°</p> <p>0 % UT (100 % dip in UT) for 250/300 cycle at 0°</p> | Mains power quality should be that of a typical commercial or hospital environment. If the user of the models TT-CL016 product name requires continued operation during power mains interruptions, it is recommended that the models TT-CL016 be powered from an uninterruptible power supply or a battery. |

| | | | |
|--|-----------------|-----------------|---|
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m, 50/60Hz | 30 A/m, 50/60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE: UT is the a. c. mains voltage prior to application of the test level. | | | |

Guidance and MANUFACTURER'S declaration - electromagnetic IMMUNITY

Guidance and manufacturer's declaration - electromagnetic immunity

The TT-CL016 is intended for use in the electromagnetic environment specified below. The customer or the user of the TT-CL016 should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level |
|-------------------------------|--|---|
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms 150 kHz to 80 MHz |
| Radiated RF IEC 61000-4-3 | 6 Vrms 150 kHz to 80 MHz outside ISM bands 10 V/m 80 MHz to 2.7 GHz | 6 Vrms 150 kHz to 80 MHz outside ISM bands 10 V/m |

Electromagnetic environment - guidance

Portable and mobile RF communications equipment should be used no closer to any part of the Models TT-CL016, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended separation distance

$$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$$

$$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$$

80MHz to 800MHz

$$d = \left[\frac{7}{E_1} \right] \sqrt{P}$$

800MHz to 2.7GHz

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these

frequency ranges.

c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TT-CL016 is used exceeds the applicable RF compliance level above, the TT-CL016 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the TT-CL016.

d. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM

Recommended separation distances between portable and mobile RF communications equipment and the model TT-CL016

The Model TT-CL016 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model TT-CL016 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model TT-CL016 as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output of transmitter W | Separation distance according to frequency of transmitter m | | |
|--|---|---|---|
| | 150 kHz to 80 MHz | 80 MHz to 800 MHz | 800 MHz to 2.7 GHz |
| | $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ | $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ | $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the

equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distances between RF wireless communications equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

| Frequency MHz | Maximum Power W | Distance | IEC 60601 Test Level | Compliance Level |
|---------------|-----------------|----------|----------------------|------------------|
| 385 | 1.8 | 0.3 | 27 | 27 |
| 450 | 2 | 0.3 | 28 | 28 |
| 710 | 0.2 | 0.3 | 9 | 9 |
| 745 | | | | |
| 780 | | | | |
| 810 | 2 | 0.3 | 28 | 28 |
| 870 | | | | |
| 930 | | | | |
| 1720 | 2 | 0.3 | 28 | 28 |
| 1845 | | | | |
| 1970 | | | | |
| 2450 | 2 | 0.3 | 28 | 28 |
| 5240 | 0.2 | 0.3 | 9 | 9 |
| 5500 | | | | |
| 5785 | | | | |

Electromag-netic Environment- Guidance

RF wireless communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended
separation distance

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARNINGS!

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transeiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

CONTRAINDICATIONS

There are no absolute contraindications to light therapy, although there are some circumstances in which caution is required. These include when the patient

- 1) has a condition that might render his or her eyes more vulnerable to phototoxicity,
- 2) has a tendency toward mania,
- 3) has a photosensitive skin condition, or
- 4) is taking a photosensitizing medication or herb (such as St. John's wort or a psoralen).

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Nota: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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FCC Compliance

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

WEEE Compliance

Correct Disposal of This Product (Waste Electrical & Electronic Equipment)
(Applicable in countries with separate collection systems)

This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take these items for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product and its electronic accessories should not be mixed with other commercial wastes for disposal.